picovest® universal NF

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Gefahr

picovest® universal NF

Phosphate-bound precision investment material for the crown and bridge technique as well as for all pressable ceramics and optimised for rapid casting

Mixing ratio: 100 g: 25 ml
150 g: 38 ml
Processing time: approx. 5 minutes

Processing

- Add liquid
- Sprinkle powder
- Mix vigorously for 20 seconds
- Stir under vacuum for 60 seconds
- Stirring speed: 250–360 rpm. —higher stirring speeds lead to less expansion
- Hold the vacuum for another 30 seconds

If investment rings are used, they must be lined with fleece as follows:

Investment ring size 1 + 3: 1 fleece insert
Investment ring size 6 + 9: 2 fleece inserts

The inside of the investment ring as well as the fleece should be greased with vaseline. picodent fleece has the same height as the investment ring. For NPM use 2 layers of fleece as of investment ring size 3, and use 3 layers of fleece as of investment ring size 6.

picovest universal NF is also suitable for muffle-free systems (e.g. picodent speed muffle system).

The casting ring is filled at the lowest vibration stage. Once the casting ring has been filled, do not vibrate again. Alternatively, place the muffle next to the vibrator and pour the investment material into the vibrator via the contact mixing cup.

Pressure pot: embedding in the pressure pot is not necessary. However, if desired, then for 10 minutes at max. 2.5 bar. Too much pressure can lead to changes in the model. Then evacuate the pressure pot slowly (approx. 4–5 minutes).

Rapid casting (recommended variant)

Preheating: after 20 minutes, calculated from the initial powder/liquid contact, the investment ring can be placed in the 850°C hot furnace.

For pressed ceramics, observe the preheating temperatures and preheating times of the respective manufacturers.

For very large resin models (e.g. picobello, picodent), the heating variant “overnight preheating” should preferably be used. For model casting with light wax, please use the picodent model casting investment material: picocast speed NF.

Preheating times

Investment ring size	1	25 minutes
	3	45 minutes
	6	75 minutes each at final temperature.

Overnight preheating

Seal the investment ring with cling film or wax lid and place in the cold furnace.

Holding times are required at 290°C for 45 minutes and at 580°C for 30 minutes.

Hold the respective final temperature for 30–45 minutes. (max. final temperature: 1050°C)

Heating speed up to 580° C at 3–5° C/min, from 580°C at 9°C/min.

Concentrate blends

50% onlays, inlays, primary components

50% gold casting, crowns and bridges (wax) and ceramic firing 60% gold casting, crowns and bridges (deep-drawing foil) and ceramic firing

70–75% secondary components (wax deep-drawing foil)

75% pressed ceramic

95–100% NPM alloy

The concentrate blends are guide values and depend on the type of alloy. For pressed ceramics – please observe the manufacturer's instructions. If expansion is too low, we recommend using the special mixing liquid for high expansion.

A mixing table with differentiated guide values for both mixing liquids is available upon request.

Storage temperature: important!! 18°–20° C are ideal.

At temperatures above 26°C, expansion can no longer be controlled correctly.

Warning: this investment material contains quartz and cristobalite. Causes damage to organs through prolonged or repeated exposure. Do not breathe dust/fume/gas/mist/vapours/spray. Get medical advice/attention if you feel unwell. Comply with local/regional/national/international regulations when disposing of contents/containers.

Danger

The above information is given to the best of our knowledge and with due diligence.

This corresponds to the current state of the art. We guarantee the flawless quality of our products, but are not liable for further processing results, which usually arise outside our sphere of influence.

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