


# DECIDE FOR NEW FLEXIBILITY

The new Implantat-rock® dental arch stone from picodent® provides you with maximum processing flexibility!

You determine the mixing ratio and the stone adapts most flexibly to the particular requirement.

Decide now for Implantat-rock®!  
It will convince you!

## Implantat-rock®



picodent® Dental-Produktions- und Vertriebs-GmbH

Lüdenscheider Str. 24-26  
51688 Wipperfurth

+49 2267 6580 0  
+49 2267 6580 30

✉ [export@picodent.de](mailto:export@picodent.de)  
🌐 [www.picodent.de](http://www.picodent.de)

## New flexibility with Implantat-rock®

Modern raw materials and the particularly fine powder structure of Implantat-rock® provide you with optimal processing conditions.

Decide the mixing ratio entirely as you wish! You have the option of variably dosing 19-20 ml of distilled water to 100 g powder and, depending on the requirement, additionally with +/- 1 ml! This way, you obtain flowability or stability simply according to your ideas.

Test the new Implantat-rock® dental arch stone now and convince yourself of its fascinating flexibility!

### FURTHER ADVANTAGES FOR YOU:

Suitable for all model systems

- ✓ Fast and high final hardness
- ✓ Super fine and creamy
- ✓ Fresh modern colours
- ✓ Low expansion
- ✓ Brilliant surface



### Technical data

Mixing ratio for	
100 g powder:	19 – 20 ml (+/- 1 ml)
Processing time:	approx. 7 min.
Stirring time under vacuum:	40 sec.
Demoulding after:	approx. 35 min.
Hardness after 24h:	approx. 300 MPa (N/mm <sup>2</sup> )
Compressive strength after 24h:	approx. 70 MPa (N/mm <sup>2</sup> )
Expansion after 2h:	< 0.08%

## ORDER NOW

### 5 KG CLIMATE BAG



1 climate bag á 5 kg  
Implantat-rock® powder

### 20 KG CARTON



4 climate bags á 5 kg  
Implantat-rock® powder

### CHOICE OF COLOURS



gold-yellow



gold-brown



ivory



grey

Discover Implantat-rock® now, the new dental arch stone with fascinating flexibility.

We look forward to you contacting us.

PHONE  
+49 2267 6580 0

FAX  
+49 2267 6580 30

INTERNET  
[www.picodent.de](http://www.picodent.de)

picodent® Dental-Produktions- und Vertriebs-GmbH

Lüdenscheider Str. 24-26  
51688 Wipperfurth

+49 2267 6580 0  
+49 2267 6580 30

[export@picodent.de](mailto:export@picodent.de)  
[www.picodent.de](http://www.picodent.de)

