pico-splint star

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
Fields of application

Fabrication of individual parts, for example, drilling templates for implant dentistry, X-ray templates, base plates, splints, fixation or transfer keys.

Important notes on processing

- **pico-splint star** can be combined with PMMA resins and deep-drawn splints (observe suitability according to manufacturer's data) for stabilising or adapting the surface.
- Observe material thicknesses, minimum 1.2 mm, maximum 5 mm, up to 10 mm for polymerisation on both sides.
- Biocompatibility is only guaranteed with complete polymerisation.
- Avoid high temperatures during polymerisation as this can lead to distortion of the objects. Non-cooled stroboscopic light units are therefore not suitable.
- Prior to polymerisation, check the contoured work for possible lifting or separating from the model.
- Close tin tightly immediately after removal of material.
- Milling and polishing should only be performed on work stations with suction devices.
- Avoid contact with mucous membranes and eyes. In case of accidental contact, rinse with adequate running water, consulting an ophthalmologist if necessary. Non-cured material can lead to irritation and sensitisation upon contact with skin or mucous membranes.
- Wear protective gloves and eye protection when processing material which has not yet set. Avoid contact with skin.
- **pico-splint star** can be disinfected using impression disinfection agents, for example.
- The product was developed for use in the dental field and may only be used by professionals for its intended use.
- Please follow the instructions on the safety data sheet.

PROCESSING TIPS:

1. Preparation

Block out plaster model, wet and carefully apply alginate isolation (repeat if required). Allow to dry thoroughly. The use of petrolatum is recommended to isolate the counter-bite in the articulator. To avoid sticking, apply a thin coat of petrolatum to fingers, gloves (i.e. nitrile), modelling instruments etc. This considerably facilitates performing and modelling.

2. Application

Remove an appropriate amount of **pico-splint star** from the tin, reclose tightly, and shape the material to a homogenous ball in your palms applying slight pressure. Assess the amount of material, homogenise and roll again after adding or removing material to avoid inclusions. Further shaping depends on the prevailing situation (i.e. denture model U-shape or basic ball). **pico-splint star** into the corresponding areas and adapt gently. Form the desired shape by gentle, repeated pressing with a finger tip. Observe appropriate dimensioning (recommended minimum thickness 1.2 mm). If the antagonists in the articulator are to be included, then apply petrolatum both to the model as well as the **pico-splint star** surface to avoid sticking and deformation. In accessible areas it is possible to smooth the surface by smoothing with finger tips.

3. Curing

pico-splint star can be polymerised in all commercially available light curing devices with a light spectrum of 320 nm - 400 nm. Units with cooling of the polymerisation chamber are also suitable. The average polymerisation time in light-curing devices with UV/UV-A lamps is approx. 10 min. Wet again for 2 - 3 before removal. This softens the insulation and facilitates removal. Perform loosening and removal of the model with great care and in direction of insertion (risk of breaking if wedged). Light-cure the opposite side of the model for approx. 5 minutes if subsequent mechanical polishing (**pico-strahresin**) is intended. If **pico-splint star Glanzlack** is used, this step is not required as a further polymerisation process is required. Please note that the polymerisation time always depends on the lamps in the device. Uniform curing times can only be ensured with regularly serviced devices. Observe the device manufacturer's instructions for use.

4. Processing

Following complete polymerisation of **pico-splint star**, remove the slightly sticky inhibition layer by rubbing with isopropanol. The fixation of functional or auxiliary parts, as well as modifications, repairs or corrections can be performed with a material approved for this application (e.g. **picopreci plus glasklar**). Please observe the instructions for use of the products used. In order to join **pico-splint star** with a deep-drawn splint, ensure that this is suitable for connecting to PMMA resins (i.e. ERKODUR, BIOCRYL). Then roughen the splint by sandblasting, apply **pico-splint star Bonder** and allow to work for 1-2 min. Then adapt **pico-splint star** and cure for approx. 10 min. Tungsten carbide burs are suited for processing and grinding. Sand the surface with emery paper and then polish mechanically. **pico-splint star Glanzlack** gloss lacquer can be used in place of manual high gloss polishing. First additionally smooth the surface by sandblasting (**pico-edelkorund 110 my**). Sealing of the top and bottom sides can be performed in two work steps. Tip: only apply the lacquer very thin to the bottom side, possibly disperse with air, so as not to impair the fit. Light-cure **pico-splint star Glanzlack** within approx. 5 min. per side.

5. Disinfection and sterilisation

Disinfection
pico-splint star can be disinfected with standard agents. Please observe the manufacturer's compatibility data.

Sterilisation

Locally applicable regulations and hygiene guidelines are to be observed. **pico-splint star** is validated for the „W&H Lisa 522, Programme UNIVERSAL 134“ (see technical specifications). Steam sterilisation may only be performed with equipment that complies with the standards EN 13060 and EN 285. Validation of the sterilisation processes was conducted according to EN ISO 17664. Responsibility for sterility lies with the user.

Technical specifications:

- **Volume/weight**
pico-splint star 60 g tin
pico-splint star Glanzlack 10 ml bottle
pico-splint star Bonder 5 ml bottle
- **Dispensing & application** Direct application, manual contouring
- **Product colour** transparent
- **Curing** Light spectrum of equipment 320 nm - 400 nm
UV/UV-A lamps approx. 10 min. on model + approx. 5 min. bottom side without model
- **Processing** at 23°C ± 2°C, 50 ± 5 % rel. humidity
- **Sterilisation process** W&H Lisa 522, Programme UNIVERSAL 134 with the following parameters:

Start	-	Heating	-
Ventilation	-0,86 bar	Sterilisation	134 -138°C; 2,10-2,20 bar; 4 min
Steam injection	+0,10 bar	Pressure relief	-
Ventilation	-0,84 bar	Drying	16 min
Steam injection	+0,50 bar	Ventilation	-
Ventilation	-0,83 bar	End	-

Contraindication

Sensitisation through **pico-splint star** is extremely rare. If such reactions occur, the application of **pico-splint star** is to be discontinued. **pico-splint star** may not be dispensed intraorally in nonpolymerised condition.

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DETAX GmbH&Co. KG
picodent Dental-Produktions- und Vertriebs GmbH