picovest[®] universal

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GB Phosphate bonded precision-quick filling investment material

Mixing ratio:	100 gr : 26 ml; 150 gr : 39 ml	Working time: approx. 5 minutes
Instructions:	• Provide requisite amount of liquid • Pour in the powder • Mix careful for 20 seconds • Mix for 90 - 120 seconds under vacuum • Stirring speed: 350-360 r.p.m – higher stirring speed leads to lower growth • Maintain the vacuum for another 30 seconds	

When using muffle rings, please use following fleece:
Muffle size 1 + 3: 1 fleece inlay
Muffle size 6 + 9: 2 fleece inlays
The muffle inside and the fleece too must be greased with vaseline. picodent fleece has the same height as the muffle.
NEM uses a 2-layer fleece per size 3 muffle, a 3-layer fleece per size 6 muffle.
picovest universal is also suitable for muffleless systems (e.g. picodent Speed Muffle System).
The padding of the casting ring is set at the lowest vibration level. Once the casting is padded, then the vibration will cease. Alternatively the muffles next to the shaker set and the dental casting investment moulds when the mixing flask gets in contact with the shaker. The casting ring will then be given up to 10 minutes max. at about 3. bar in the accumulator, so that the accrued pressure may lead to changes in the temporary cast. Then slowly empty the accumulator (ca. 4-5 minutes).
Quick casting
Preheating: 20-25 minutes after the mixing, the casting ring can be set in the oven heated up at 850°C. Observe each manufacturer's preheating temperatures and preheating times for press ceramics. For larger acrylic casts, it is preferable to set the heating option to „Overnight preheating“. For castings with light wax, please use the Picodent dental casting investment.
Preheating times
Muffle size 1x25 minutes 3x45 minutes 6x75 minutes meant for end temperature
Over night heating
Seal the former with wax or fresh keeping foil and set in cold oven. Waiting time is required on 290°C for 45 minutes and on 580°C for 30 minutes.
Keep the end temperature for 30-45 minutes. (maximum end temperature: 1050°C). The speed of heating until 580°C in 3-5°C/minute, from 580°C in 9°C/minute.

Concentrate mixes
50% onlays, inlays, primary parts
50% cast gold, crowns and bridges (wax) and baked-on ceramics
60% cast gold, crowns and bridges (thermoforming foil) and baked-on ceramics
70-75% secondary parts (wax-thermoforming foil)
75% press ceramics
95-100% NEM alloy
The concentrate mixes are guideline values and depend on the type of alloy. In the case of press ceramics – please follow the manufacturer's instructions. If expansion rate is insufficient, we recommend the use of the special mixing liquid for high expansion.
A mixing table with differentiated benchmarks can be requested.
Storage temperature: Important!! Storage at 18° – 20°C is ideal. For temperatures over 26°C, the expansion can no longer be managed with accuracy.
Warning: these investment material contains quartz and cristobalite. May cause damage to the lung through prolonged or repeated exposure. Route of exposure: Inhalative. Do not breathe dust. In case of inadequate ventilation wear respiratory protection. Get medical advice/attention if you feel unwell.
The recommendations are given to the best of our knowledge after careful control. We guarantee the quality of our products. Any further liability cannot be accepted since the proper application of our products is outside of our control.

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