twinduo scan

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Farbe: gelb
Schraubdeckeldose: 2 x 1,5 kg

twinduo scan

GB Scanable two-component silicone putty, especially designed for use in CAD/CAM systems

for dentistry

Shore hardness (10 min): 90
End hardness: 95

User manual

Dose: 1 measuring spoon of lab putty (base), yellow - 1 measuring spoon of lab putty (catalyst), grey

Preparation: 1 spoon full of yellow base material (approx. 8 ml/14 g) is mixed with a spoon full of grey catalyst material (approx. 8 ml/14 g) using the hands. The different colour of the base and catalyst materials means that it is possible to see when the substances have been well mixed.

Mixing time: 30 seconds
Total working time: 1.5 minutes
Setting time: 5-6 minutes

These values are based on a working temperature of 23°C

The hardened twinduo is boil-proof and maintains its forms when used in the microwave (600 Watt, 850°C).

Attention: The material heats up to 60°C after just 10 seconds in the microwave.

Low temperatures and lower doses of the catalyst lengthen the setting time, and high temperatures and higher doses of the catalyst shorten the setting time.

Use by: see packaging label, at least 2 years

Storage: Do not store over 25°C.

Working and safety instructions:

- When mixing by hand **do not use latex** gloves
- Close the **box** immediately after use
- **Do not mix up the dosing spoon and the lids**
- the castings can be copper-plated or silver-plated using galvanisation
- Avoid contact with the eyes, but if the material comes into contact with the eyes, rinse them immediately with plenty of clean water, and seek medical assistance straight away. Avoid splashing the clothing, since curing elastomer setting materials are resistant to chemicals.

Delivery forms: Colour: yellow - Box with screw-on lid: 2 x 1.5 kg

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