picocast speed NF

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picocast speed NF

Investment material for the model casting technique and optimised for speed-casting

Mixing ratio:
For silicone: 100 g powder : 20 ml liquid 200 g powder : 40 ml liquid 400 g powder: 80 ml liquid

- Processing**
- Add liquid.
 - Sprinkle powder.
 - Mix vigorously for 20 seconds.
 - Stir under vacuum for 60 seconds.
 - Stirring rate: 250-360 rpm
 - Hold the vacuum for another 15-20 seconds.

The casting ring is filled at the lowest vibration level or the ring is placed next to the vibrator and the investment material is poured in with the vibrator. Once the casting ring has been filled, do not vibrate again. Remove the muffle former and roughen it on the opposite side to the funnel (make three scratches).

Processing time: approx. 5 min.
Setting end: approx. 9 min.

To achieve optimum casting results, we recommend a 50% liquid concentration for overinvestment.
Pressure pot: embedding in the pressure pot should be up to max. 2.5 bar. It can remain there until max. before mounting.

Liquid concentration for the model

75-80 % normal model casting
80-85 % light waxes
75-80 % combination work, depending on the size and extent of the jaw
The concentrate blends are guide values and depend on the type of alloy.

Special recommendation for light waxes

450 g powder (9er investment ring): 59 ml liquid / 31 ml dist. water
400 g powder: 52 ml liquid / 28 ml dist. water
300 g powder (6er investment ring): 39 ml liquid / 21 ml dist. water

Speed-casting technique
(recommended variant)

25 minutes after filling, the investment ring can be placed in the 850-950°C hot furnace.
For furnaces with floor heating, keep sufficient distance (approx. 1 cm) between investment ring and floor plate.
Maintain the alloy-dependent final temperature for at least 60 minutes up to 400 g investment ring weight and for at least 90 minutes up to 600 g investment ring weight.

Overnight pre-heating

Place the muffle in the cold furnace after curing. At 280°C and at 580°C, a holding time of 45-60 minutes is required, depending on the size and number of investment rings.
Maintain the alloy-dependent final temperature for at least 60 minutes up to 400 g investment ring weight and for at least 90 minutes up to 600 g investment ring weight.

Heating rate approx. 3-7°C/min. at linear furnace control.

Maximum final temperature 1030°C. Pour rapidly, avoid casting delays!

Warning
This investment material contains quartz and cristobalite.
Causes damage to organs through prolonged or repeated exposure.
Do not breathe dust/fume/gas/mist/vapours/spray. Get medical advice/attention if you feel unwell.
Dispose of contents / container in accordance with local/regional/national/international regulations.

Do not open the pre-heating furnace during the heating phase, as the wax vapours produced can ignite in the air!
The above information is given to the best of our knowledge and with due diligence. This corresponds to the current state of the art. We guarantee the flawless quality of our products, but are not liable for further processing results, which usually arise outside our sphere of influence.

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