# picovest<sup>®</sup> royal NF

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Gefahr

# picovest<sup>®</sup> royal NF

Graphite-free phosphate-bound investment material for the crown and bridge technique, especially suited for NPM and optimised for rapid casting

**Mixing ratio:** 100 g : 24 ml  
150 g : 36 ml

**Processing time:** approx. 5 minutes

#### Processing

- Add liquid
- Sprinkle powder
- Mix thoroughly with a spatula by hand for 30 seconds
- Stir under vacuum for 60 seconds / stirring rate 250–350 rpm

It is preferable to use **picovest<sup>®</sup> royal NF** with the ringless **picodent<sup>®</sup> Speed muffs**. Metal muffle rings can also be used. Fill the prepared muffle former evenly with the investment material. If the object has a complicated geometry, the investment material can be poured in under slight vibration (vibrator at the lowest level). Excessive vibration must be avoided as this may cause bubble formation and segregation.

**The investment material does not need to be cured under pressure.**

#### Rapid casting (recommended variant)

- Fill the muffle former with investment material and let it set for 15 minutes (time starts from initial contact powder/liquid).
- Remove the muffle former and roughen it on the opposite side to the funnel (make three scratches).
- Place the muffle in the furnace preheated to a maximum of 850 °C no later than 15 minutes after powder-liquid contact. Holding times at 850 °C: Size 2: 30 min. / Size 3: 45 min. / Size 6: 60 min. / Size 9: 90 min.
- Depending on the alloy used, heat the furnace to the final temperature (950 °C max.). In any case, the holding times at 850 °C should be adhered to.

Always apply a wax coat to resin models (modelling resin, self-curing/powder-liquid or **picobello plus**/light-curing).

#### Slow heating (conventional casting)

- Fill the muffle former with investment material and let it set for 20 minutes.
- Remove the muffle former and roughen it on the opposite side to the funnel (make three scratches).
- Place the muffle in the cold furnace. For larger resin models, program a holding stage at 250 °C. Size 2: 30 min. / Size 3: 45 min. / Size 6: 60 min. / Size 9: 90 min
- Heat the furnace evenly to the final temperature of 850°–max. 950°C. For this also, please follow the instructions for the alloy being used. Recommended temperature increase rate: 9°C/min. A further holding stage is not required.
- Let the muffle stand in the furnace at the final temperature depending on the muffle size used (size 2: 30 min. / size 3: 45 min. / size 6: 60 min. / size 9: 90 min.).
- Tip: With overnight preheating, seal the muffle with a cling wrap or wax lid.

#### Concentration blends for 150g powder

Indication	Modelling material	Casting alloy	Concentration [ml] liquid: aqua dest	Concentration [%]	In-house values
Crowns & bridges	Wax	Nonprecious metal alloys	32 : 4	90	
		Reduced-gold (> 55%) and metal-ceramic alloys for porcelain	25 : 11	70	
		Palladium-based alloys	27 : 9	75	
		High-gold content (> 70% Au)	23 : 13	65	
Inlay (trihedral)	Wax	Nonprecious metal alloys	32 : 4	90	
		Reduced-gold (> 55%) and metal-ceramic alloys for porcelain	25 : 11	70	
		Palladium-based alloys	27 : 9	75	
		High-gold content (> 70% Au)	23 : 13	65	
Female parts, telescopic and conical crowns 0-6°	Modelling resin, self-curing (powder/liquid)	Nonprecious metal alloys	32 : 4	90	
		Reduced-gold (> 55%) and metal-ceramic alloys for porcelain	22 : 14	60	
		Palladium-based alloys	22 : 14	60	
		Hochgoldhaltig (> 70% Au)	20 : 16	55	
Female parts, telescopic and conical crowns 0-6°	picobello plus, light-curing with wax coat	Nonprecious metal alloys	34 : 2	95	
		Reduced-gold (> 55%) and metal-ceramic alloys for porcelain	23 : 13	65	
		Palladium-based alloys	22 : 14	60	
		High-gold content (> 70% Au)	22 : 14	60	

The concentrations given are merely for guidance. These can be adjusted by ± 5 %.

All information is based on our laboratory's test results and considered a guideline. All information is subject to change. Please also follow the directions for use of the alloy used.

**Storage and processing temperature: 20-23°C is ideal!**

**As a general rule of thumb, consistent storage and processing = consistent results**

**Warning:** this investment material contains quartz and cristobalite. Causes damage to organs through prolonged or repeated exposure. Do not breathe dust/fume/gas/mist/vapours/spray. Get medical advice/attention if you feel unwell. Comply with local/regional/ national/ international regulations when disposing of contents/containers.

The above information is given to the best of our knowledge and with due diligence. This corresponds to the current state of the art. We guarantee the flawless quality of our products, but are not liable for further processing results, which usually arise outside our sphere of influence.

Please also observe the relevant valid safety data sheets!

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