IT

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picopreci plus trasparente

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Light-Curing Precision Resin

Description

picopreci plus transparent is a light-curing gel-like precision modelling resin. It makes precisely accurate work possible both when used on models in a dental laboratory and directly in a patient's mouth. **picopreci plus transparent** polymerises in all standard light-curing units with a light spectrum of 360-420 nm. Curing is done without vacuum.

Composition

Preparation of acrylates and phosphine oxides.

Indication

- Repairing, building up and correcting adjusted bite guard splints
- Positioning and fixing of functional parts to orthodontic appliances, radiographic templates and drilling templates
- Fabrication of fixation and implant jigs
- Fine corrections in resin (e.g. papillae)

Characteristics and Benefits

picopreci plus makes fast and precise work possible with short polymerisation times. The material is characterised by precise application directly from the syringe with a dispensing tip and by its high stability, i.e. it requires little use or time, and placement is exact.

picopreci plus transparent is methyl methacrylate-free.

Instructions for Use

After removing the sealing cap, you should attach the enclosed dispensing tips for precisely accurate application. These should only be used once. Reseal the push syringe afterwards. Aim for a consistent layering thickness (about 1-2 mm) during the application. Light units with stroboscopic light (e.g. Dentacolor) require 3 minutes, while light units with UV-A tubes (e.g. Megalight) need 10-20 minutes depending on the number of tubes. Bright plaster models (Implant-rock, picodent) assist curing.

Curing with LED lamps is possible in 20 seconds if their emission range is 360-420 nm (e.g. pico-light LED).

Please note that the curing time is critically dependent on the emission performance of the lamps used. Correct maintenance of the devices and regular replacement of the light fixtures contribute to consistent curing. Also follow the particular device manufacturer's instructions.

If using in a patient's mouth, please be aware that the material heats up during curing. To avoid unacceptable temperature peaks, you should therefore build up larger amounts (> 0.5 g) in layers and polymerise.

Important: When working on the model, always wet with alginate insulation agent (plaster/plastic insulation agent, picodent). With application to PMMA plastics for chemical retention, always dab the corresponding areas outside the mouth with **picopreci plus Bonder** beforehand and light-cure for about 3 minutes. For mechanical retention, the parts can additionally be prepared by being roughened with sandpaper or blasting.

To be able to apply material after curing, do not remove the inhibition layer on the surface; it ensures optimum bonding and prevents stresses in the modelling resin. If the inhibition layer has already been removed, **picopreci plus Bonder** should be applied again for reliable bonding.

Only after final completion is the surface released from the inhibition layer and polished with rotating instruments or isopropyl.

Warnings

Please also observe the valid safety data sheets.

picopreci plus transparent contains TPO, UDMA.

May cause an allergic skin reaction. Toxic to aquatic life with long lasting effects. Avoid breathing dust/fume/gas/mist/vapours/spray. Wear protective gloves/protective clothing/eye protection/face protection.

IF ON SKIN: Gently wash with plenty of soap and water. If skin irritation or rash occurs: Get medical advice/attention.

Bonder contains: MMA, PETIA/PETTA, TPO. Highly flammable liquid and vapour. Causes skin irritation. May cause an allergic skin reaction. Causes serious eye damage. May cause respiratory irritation. Harmful to aquatic life with long lasting effects. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Avoid breathing dust/fume/gas/mist/vapours/spray. Avoid release to the environment. Wear protective gloves/protective clothing/eye protection/face protection. IF ON SKIN: Gently wash with plenty of soap and water.

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Switch on an exhaust fan during trimming.

Contraindications, Interactions

Do not use the product if you have a known allergy to one of its materials. Otherwise patient group not restricted. There are no known Interactions.

Secondary Effects

Product may cause allergic reaction.

Storage Instructions

Store at 4-25 °C. Keep **picopreci plus** syringes away from heat and sunlight. Do not use after the expiry date.


picopreci plus Set: 3 x 3 g **picopreci plus** colour: **transparent** in light-protected Slide syringes incl. 12 dispensing tips, 5 ml **picopreci plus Bonder**

picopreci plus Set: 5 x 3 g **picopreci plus** colour: **transparent** in light-protected Slide syringes incl. 20 dispensing tips

The product was developed for use in the dental field and must be used in accordance with the directions for use. The manufacturer accepts no liability for damage caused by the product being used in other ways. Furthermore, the user is obliged to test the material on his or her own responsibility for suitability and usability for the intended purposes before it is used, especially if this is not specified in the directions for use. Serious incidents must be reported to DeltaMed GmbH and the responsible authority.

Disposal Instructions

Dispose of contents/container in accordance with local regulation.

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picodent Dental-Produktions- und Vertriebs GmbH