picovest® royal

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picovest® royal

Graphite-free phosphate-bound investment for the crown and bridge technique for conventional and fast heating

Mixing ratio: 150 g: 36 ml
Processing time: approx. 6 minutes

- Processing**
- Add liquid
 - Sprinkle powder
 - Mix thoroughly with a spatula by hand for 30 seconds
 - Stir under vacuum for 60 seconds / stirring rate 350 RPM

It is preferable to use **picovest® royal** with the ringless **picodent® Speed muffles**. Metal muffle rings can also be used.

Fill the prepared muffle former evenly with the investment. If the object has a complicated geometry, the investment can be poured in under slight vibration (vibrator at the lowest level). Excessive vibration must be avoided as this may cause bubble formation and segregation.

The investment does not need to be cured under pressure.

- Slow heating (conventional casting)**
- Fill the muffle former with investment and let it set for 20 minutes.
 - Remove the muffle former and roughen it on the opposite side to the funnel (make three scratches).
 - Place the muffle in a cold oven. For larger resin carvings, program a retention shoulder at 250 °C. Size 2: 30 min. / Size 3: 45 min. / Size 6: 60 min. / Size 9: 90 min.
 - Heat the oven evenly to the final temperature 850° - max. 950°C. For this also, please follow the instructions for the alloy being used. Recommended temperature increase rate: 9°C/min. **A further holding stage is not required.**
 - Let the muffle stand in the oven at the final temperature depending on the muffle size used (size 2: 30 min. / size 3: 45 min. / size 6: 60 min. / size 9: 90 min.).

→ Tip: With **overnight preheating**, seal the muffle with a cling wrap or wax lid.

- Fast heating (speed casting)**
- Fill the muffle former with investment and let it set for 20 minutes.
 - Remove the muffle former and roughen it on the opposite side to the funnel (make three scratches).
 - Place the muffle in an oven preheated to a maximum of 850 °C no later than 25 minutes after powder-liquid contact. Holding times at 850 °C: Size 2: 30 min. / Size 3: 45 min. / Size 6: 60 min. / Size 9: 90 min., **always position resin carvings at 600 °C.**
 - Depending on the alloy used, heat the oven to the final temperature (950 °C max.). In any case, the holding times at 850°C should be adhered to.
 - Let the muffle stand in the oven at the final temperature depending on the muffle size used (size 2: 30 min. / size 3: 45 min. / size 6: 60 min. / size 9: 90 min.).

Always apply a wax coat to resin carvings (modelling resin, self-curing/powder-liquid or **picobello plus**/light-curing).

Concentration Recommendations for 150 g Powder

Indication	Carving material	Casting alloy	Concentration [ml] Liquid: aqua-dest	Concentration [%] Liquid: aqua-dest	In-house values
Crowns & bridges	Wax	Nonprecious metal alloys	30 : 6	83 : 17	
		Reduced-gold (> 55%) and metal-ceramic alloys for porcelain	21 : 15	58 : 42	
		Palladium-based alloys	23 : 13	64 : 36	
		High-gold (> 70% Au)	20 : 16	56 : 44	
Inlay (trihebral)	Wax	Nonprecious metal alloys	29 : 7	81 : 19	
		Reduced-gold (> 55%) and metal-ceramic alloys for porcelain	21 : 15	58 : 42	
		Palladium-based alloys	23 : 13	64 : 36	
		High-gold (> 70% Au)	20 : 16	56 : 44	
Female parts, telescopic* and conical crowns 0-6°	Modelling resin, self-curing (powder/liquid) with wax coat	Nonprecious metal alloys	32 : 4	89 : 11	
		Reduced-gold (> 55%) and metal-ceramic alloys for porcelain	17 : 19	47 : 53	
		Palladium-based alloys	20 : 16	56 : 44	
		High-gold (> 70% Au)	15 : 21	42 : 58	
Female parts, telescopic* and conical crowns 0-6°	picobello plus, light-curing with wax coat	Nonprecious metal alloys	33 : 3	92 : 8	
		Reduced-gold (> 55%) and metal-ceramic alloys for porcelain	18 : 18	50 : 50	
		Palladium-based alloys	21 : 15	58 : 42	
		High-gold (> 70% Au)	16 : 20	44 : 56	

* With telescope crowns, the concentration can be increased by +1 ml liquid / -1 ml distilled water (36 ml overall must be adhered to).

All information is based on our laboratory's test results and considered a guideline. All information is subject to change. Please also follow the directions for use of the alloy used. Individual mixing concentrations can be entered in the empty table.

Storage and processing temperature: 20-23°C is ideal. As a general rule of thumb, consistent storage and processing = consistent results!

Warnings

Warning

picovest® royal contains quartz and cristobalite. Appropriate suction or ventilation must be ensured at the processing site. Inhaling dust from the product may cause silicosis. P2 dust masks must be provided by the employer. Ammonia is released when the investment is heated. Ammonia is a respiratory irritant. Skin and mucous membrane irritations (mechanical) are possible in the respiratory tract and eyes. Localised contact is likely to cause irritation and potentially corrosivity to the skin and mucous membranes.

picovest® royal Liquid: If contact with the skin occurs: Rinse with plenty of water. If contact with the eyes occurs: Hold the eye open and rinse with plenty of clean water. If symptoms persist, consult an ophthalmologist. Upon swallowing: Rinse the mouth with water and drink plenty of water. Never pour anything into the mouth of an unconscious person. Mixing liquids are alkaline.

Please also observe the relevant valid safety data sheets.

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