dentabs

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GB The pill for gypsum separators

trimming does not settle in conventional gypsum collector systems and leads to pipe blockage. It is bad enough if a pipe branch is blocked. However, things become extreme when collecting pipes are affected and nothing in the whole building works anymore.

Prevent pipe blockages with dentabs

If a sewage pipe gets blocked, then you have a problem. And the costs add up. Downtime throughout the entire laboratory, delayed delivery, annoyed customers or even losing customers, in conjunction with flooding, a bad smell and unplanned costs, this all makes pipe cleaning anything but pleasant and inexpensive.

And this is how dentabs works

The **dentabs** tablet was specially developed for use in sedimentation cleaners in dental laboratories. Its active components bind the suspended matter from gypsum, quartz and phosphate-bonded investment materials, thus enabling complete separation in the separators. The clogging of waste water pipes by solid particles is eliminated if the separator system is emptied on time.

After emptying, simply hang the **dentabs** in the plaster separator with the supplied dosing cage (in the case of multi-chamber separators preferably at the overflow to the second chamber) – so that the tablet is completely below the water surface and is not directly exposed to the gypsum-loaded water inflow. **dentabs** soon begins to release its active ingredients. Even the finest suspended gypsum solids are flocculated and sink to the ground, leaving only clear water to drain off.

dentabs is environment-friendly

DENTABS are neutral and harmless to use. The active ingredients of DENTABS naturally comply with the applicable wastewater management regulations.

Separator cleaning made easy

Experience it yourself: emptying has never been easier! **dentabs** prevents plate formation. And where nothing hardens, there is also nothing to scrape.

Hygienic and odourless

Anyone with a nose for problems will appreciate this – because clean water does not smell and suspended solids are bound. But discipline is still required: beverage residues and soapy suds do not belong in your separator.

How often do you need dentabs?

Quite simple: one pill is enough for the gypsum capacity of a standard separator (approx. 15 litres). Best to empty regularly once a week, which means, one **dentabs** every week. You save time, you save money, and you are doing something for your employees and the environment. For yourself and for your wallet.

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