pico-print model2

385 nm / 405 nm

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
pico-print model2

385 nm / 405 nm

Light-curing material for the fabrication of dental working models.

Important note regarding use

Please read these instructions carefully and in their entirety before using the product for the first time. The product must only be used as described in these instructions. Do not use the product for any purpose other than that specified. Please follow all safety precautions and instructions contained herein. The product must only be used by persons who have read and understood these instructions.

Processing

- During the application of pico-print model2 we recommend wearing personal protective equipment e. g. suitable gloves, safety goggles etc.
- pico-print model2 should be well shaken approx. 1 hour before use.
- Carefully pour pico-print model2 into the pre-determined container of the production unit.
- Please see further steps within printer manual.
- Remove all bubbles with a cleaned object.

Processing temperature

Sonic XL 4K (qualified by Dreve) /	
Sonic XL 4K 2022 (qualified by Dreve)	23 ± 3 °C
Rapid Shape DII Drucker	23 ± 3 °C
Asiga MAX UV	35 ± 3 °C

- Select the product specific parameter settings for pico-print model2 for your process.
- When the building process is finished a direct post treatment is recommended.
- After the platform is taken up a dripping off time of approx. 10 minutes is recommended.
- Resin coated parts get clean with isopropanol (97 %) within a time of about 2 x 6 minutes. For a particular effective cleaning with isopropanol the use of explosion-protected ultrasonic units is recommended. Alternatively, cleaning with isopropanol can be carried out in a separate beaker under pouring conditions for approx. 2 x 6 minutes. Recommendation: Beaker and magnetic stirrer, laboratory shaker or ultrasonic cleaner. In case of insufficient cleaning (such as shiny areas on the object after drying) repeat cleaning process.

Post curing pico-print model2

beige-opaque

PCU LED N ₂	
Rapid Shape DII /	
Asiga MAX UV	
Sonic XL 4K (2022)	8 minutes, 80 % light power (vacuum)
FotoDent® flash	
Rapid Shape DII /	8 minutes + 1 minute gas pre-flow /
Asiga MAX UV	12 minutes + 1 minute gas pre-flow /
PCU 90	8 minutes (incl. 4 minutes gas pre-flow)
RS Cure	According to respective parameter set for pico-print model2

- Impurity due to the pico-print model2. It is recommended to regularly take out the container of the production unit, to homogenize and filtrate the content. Upcoming entrained bubbles get eliminated through a holding time of approx. 1–2 hour with pico-print model2.


To avoid damage of material quality do not expose the liquid material to light under any circumstances. Variations from the described manufacturing process may lead to modified mechanical characteristics and/or colour variations of the pico-print model2 material.

Hazards identification

Hazard and safety notes can be found on the product label and the corresponding safety data sheet.

Lot number / Expiry date

The lot number and the expiry date are indicated on each pico-print model2 packaging. In case of claims please always indicate the lot number of the product. Do not use the product after the expiry date.

Rev. 2026-06