



DaVinci's *Majesthetik*® - Stumpfix®

Dual-component material

Majesthetik®-Stumpfix® is a unique dual-curing material. After pick-up impressions, precise resin dies can be fabricated from zirconia or metals frames (precious/non-precious).

Due to the materials' light curing properties, Majesthetik® retention pins can be fixed in the lamp (Majesthetik® Fixierungslampe). Afterwards, the resin cures by itself and holding the pins is not necessary anymore. Majesthetik®-Stumpfix® is a modern system for easy and rapid processing.



Get more information and
a short video of application.
www.picodent.de/stumpfix


Item no.	Item	Contents
1550 4220	Majesthetik®-Stumpfix® Starter-Set	2 x 16 g dual-chamber syringes, 16 x Majesthetik® Automix mixing tips, 1 x Majesthetik® Fixierungslampe, Separator 50 ml + Brush, 100 x Majesthetik® Retentionspins, 1 Majesthetik®-Stumpfix® Plexiglasständer
1550 4216	Majesthetik®-Stumpfix® Set	2 x 16 g dual-chamber syringes, 16 x Majesthetik® Automix mixing tips
1550 4516	Majesthetik®-Stumpfix® Refill pack	5 x 16 g dual-chamber syringes
1550 4550	Majesthetik®-Stumpfix® Separator	50 ml + Brush
1550 4040	Majesthetik®-Stumpfix® Automix-Kanülen	Refill pack: 40 Automix mixing tips
1550 0300	Majesthetik® Fixierungslampe	1 piece, polymerisation lamp
1550 4520	Majesthetik®-Stumpfix® Plexiglasständer	1 piece plexiglas rack

DaVinci's *Majesthetik*® - Retentionspin

For resin dies after pick-up impressions

Guarantees perfect retention to dies and models. Particularly well suited for thin dies.

Item no.	Item	Contents
1550 1100	Majesthetik®-Retentionspin	100 pieces
1550 1500	Majesthetik®-Retentionspin	500 pieces



DaVinci's **Majesthetik® - Gingi**implant

Gingival mask specifically for implant technique
 With Majesthetik®-Gingiimplant, the user has the possibility to reconstruct the gum line.
 Extreme hardness (shore hardness: 70). Can be grinded very well, natural colours, high stability.

Item no.	Item	Contents
1550 2250	Majesthetik®-Gingiimplant	Set with 2 x 50 ml double cartridge including 12 mixing tips
1550 2850	Majesthetik®-Gingiimplant	Refill pack: 8 x 50 ml double cartridge
1550 5040	Majesthetik®-Zahnfleischmaske Mischkanülen	Refill pack: 40 mixing tips

DaVinci's **Majesthetik® - Gingiva**

Gingival mask for the marginal areas of crowns and bridges
 Majesthetik®-Gingiva combines perfect properties: opaque colour, high stability, high tensile strength and high final hardness (shore hardness: 62).

Item no.	Item	Contents
1550 1250	Majesthetik®-Gingiva	Set with 2 x 50 ml double cartridge including 12 mixing tips
1550 1850	Majesthetik®-Gingiva	Refill pack: 8 x 50 ml double cartridge
1550 5040	Majesthetik®-Zahnfleischmaske Mischkanülen	Refill pack: 40 mixing tips

DaVinci's **Majesthetik® - Gingiva**scan

Gingival mask specifically designed for use in CAD/CAM technology
 Scan version for powder-free, digital 3D data capture.
 High final hardness (Shore hardness: 70) and very good stability.

Item no.	Item	Contents
1550 3250	Majesthetik®-Gingiva scan	Set with 2 x 50 ml double cartridge including 12 mixing tips
1550 5040	Majesthetik®-Zahnfleischmaske Mischkanülen	Refill pack: 40 mixing tips

Separator for Gingival mask see page 27: Majesthetik-Separator for spraying and brushing.



DaVinci's **Majesthetik®-Separator**

Separating agent for spraying or brushing

Separating agent perfectly suited for Majesthetik® gingival masks against A or C silicones and for separating polyether castings before casting with polyurethane (e.g. picopoly).

Item no.	Item	Contents
1550 0050	Majesthetik®-Separator	50 ml bottle with brush top and vaporizer



DaVinci's **Majesthetik®-Line**

picodent® injector for double cartridges DS 50

For the application of Majesthetik®-Gingiva, Majesthetik®-Gingiimplant and esthetic-gum in double cartridges DS 50.

Item no.	Item	Contents
1550 1002	picodent® Injektor	1 piece



DaVinci's **Majesthetik®-Texturpuder**

To help distinguish anatomical structures

Outlines and surface textures of the teeth become more perceptible. Can be steamed off easily.

The pasty texture powder is placed on the crown surface with a brush. Simple handling - easy to remove.

Orange fragrance

Item no.	Item	Contents
1550 0200	Majesthetik®-Texturpuder silver, incl. brush	20 g glas
1550 0250	Majesthetik®-Texturpuder gold, incl. brush	20 g glas



DaVinci's *Majesthetik®-Laseal*

Light-curing stone glaze
To harden and seal stone dies.

DaVinci's *Majesthetik®-Divi-Fluid K*

Separator for stone / ceramic
Special separator for ceramics.
Very economical use, good separating effect.

Item no.	Item	Contents
1550 0010	Majesthetik®-Laseal	10 ml bottle with brush top

Item no.	Item	Contents
1550 1012	Majesthetik®-Divi-Fluid K	12 ml bottle with brush top