



## pico-edelkorund

Fine blasting medium for model casting and ceramic alloys

aluminum oxide:

- minimum of 99.8 % aluminum oxide
- free from iron and quartz
- high final hardness, remains sharp-edged (Investment material clears away quickly)
- not hazardous to health (no risk of silicosis)

Item no.	Item	Contents
1250 5005	pico-edelkorund, 50 $\mu$	5 kg dosage bag
1250 5020	pico-edelkorund, 50 $\mu$	4 x 5 kg (20 kg) dosage bag
1250 1105	pico-edelkorund, 110 $\mu$	5 kg dosage bag
1250 1120	pico-edelkorund, 110 $\mu$	4 x 5 kg (20 kg) dosage bag
1250 2505	pico-edelkorund, 250 $\mu$	5 kg dosage bag
1250 2520	pico-edelkorund, 250 $\mu$	4 x 5 kg (20 kg) dosage bag



## pico-glanzstrahlperlen

For cleaning and smoothing of precious metal alloys

glass polishing beads

- silicon oxide, super fine
- produces smooth, silky finish
- gentle on delicate surfaces
- condenses metal surfaces

Item no.	Item	Contents
1250 0505	pico-glanzstrahlperlen, 50 $\mu$	5 kg dosage bag
1250 0520	pico-glanzstrahlperlen, 50 $\mu$	4 x 5 kg (20 kg) dosage bag
1250 1205	pico-glanzstrahlperlen, 125 $\mu$	5 kg dosage bag
1250 1220	pico-glanzstrahlperlen, 125 $\mu$	4 x 5 kg (20 kg) dosage bag



## pico-bimssteinpulver

For polishing and treatment of acrylic dentures

pumice powder

- free from quartz and germs
- skin-friendly, with neutral smell
- not hazardous to health (no risk of silicosis)

Item no.	Item	Contents
1250 0205	pico-bimssteinpulver, feine Körnung (fine grain size)	5 kg dosage bag
1250 0220	pico-bimssteinpulver, feine Körnung (fine grain size)	4 x 5 kg (20 kg) dosage bag
1250 0005	pico-bimssteinpulver, mittlere Körnung (medium grain size)	5 kg dosage bag
1250 0020	pico-bimssteinpulver, mittlere Körnung (medium grain size)	4 x 5 kg (20 kg) dosage bag
1250 0105	pico-bimssteinpulver, grobe Körnung (coarse grain size)	5 kg dosage bag
1250 0120	pico-bimssteinpulver, grobe Körnung (coarse grain size)	4 x 5 kg (20 kg) dosage bag



## pico-strahlresin

Designed especially for polishing acrylic dentures

Polishing agent especially for acrylics.

Optimum polishing results - conventional processing.

- unique crystalline structure
- gentle, but extremely abrasive
- smoothing effect and longer service life (30%)
- prevents bacterial growth and is odourless

Item no.	Item	Contents
1250 3002	pico-strahlresin	2 kg dosage bag
1250 3008	pico-strahlresin	4 x 2 kg dosage bag
1250 3010	pico-strahlresin	10 kg box



## pico-bims-clean

Ready-to-use mixing liquid for pumice powder and pico-strahlresin

- high bactericidal and fungicidal effect
- effectively prevents the development of unpleasant odours
- continuously delays dehydration of the polishing agent
- excellent material compatibility

Item no.	Item	Contents
1250 4001	pico-bims-clean	1.000 ml bottle
1250 4005	pico-bims-clean	5.000 ml canister