

eco-sil

Addition-curing duplicating silicone



Addition-curing duplicating silicone for dentistry, mixing ratio 1:1

Shore A hardness: 22 / Working time: 6-7 minutes
Demouldable: after 25 min. / Hardened colour: lilac



YOUR BENEFITS AT A GLANCE:

- ✓ Very low viscosity
- ✓ High tearing strength
- ✓ Extremely sharp details due to lilac colour
- ✓ Suitable for all duplications

AVAILABLE AS:

- 2 x 1 kg bottle per component A and B (2 kg)
- 2 x 5 kg canister per component A and B (10 kg)

eco-sil speed

Addition-curing duplicating silicone - fast setting



Addition-curing duplicating silicone for dentistry, mixing ratio 1:1

Shore A hardness: 22 / Working time: 3-4 minutes
Demouldable: after 12 min. / Hardened colour: lilac



YOUR BENEFITS AT A GLANCE:

- ✓ Enables very fast working
- ✓ Can already be demoulded after 12 minutes
- ✓ Very low viscosity
- ✓ High tearing strength
- ✓ Extremely sharp details due to lilac colour
- ✓ Suitable for all duplications

AVAILABLE AS:

- 2 x 1 kg bottle per component A and B (2 kg)
- 2 x 5 kg canister per component A and B (10 kg)