



### pico-splint star

Light-curing, transparent modelling resin

Fabrication of individual parts, for example, drilling templates for implant dentistry, X-ray templates, base plates, splints, fixation or transfer keys.

Item no.	Item	Contents
1820 2106	pico-splint star Paste	60 g tin
1820 2105	pico-splint star Primer	5 ml bottle
1820 2110	pico-splint star Lack	10 ml bottle
1820 2120	pico-splint star Einmalpinsel	100 pieces blue metallic



### picobello plus gel, crystal blue

Light-curing modelling resin

A new generation of light-curing modelling resin for dental technology applications, especially for milling and attachment work.

Patterns can be tried out in the mouth prior to moulding or milling. The high translucency of picobello plus gel allows the thickness of the layers to be checked more easily, especially when modelling and finalising.

picobello plus gel can be applied directly from the dosage syringe.

Item no.	Item	Contents
1820 0040	picobello plus gel Set, crystal blue	3 x 3 g picobello plus gel in light-protected syringes, incl. 12 application tips
1820 0050	picobello plus gel refill pack, crystal blue	5 x 3 g picobello plus gel in light-protected syringes, incl. 20 application tips
1809 0012	picobello plus gel Dosierspitzen	refill pack: 50 application tips



### picobello gel

Light-curing acrylic resin for modelling

Light-curing modelling resingel supplied in dosing syringe, particularly for milling and attachment work, crowns and bridges and the interlocking of cast objects on the model and in the mouth, etc. – flowing application with high stability. picobello gel und picobello form are compatible and can therefore be used very well together.

Item no.	Item	Contents
1820 0010	picobello gel Set	3 x 3 g picobello in light-protected syringes incl. 12 application tips
1820 0020	picobello gel Set	refill pack: 5 x 3 g picobello in light-protected syringes incl. 20 application tips
1809 0012	picobello Dosierspitzen	refill pack: 50 application tips



## picopreci plus

### Light-curing precision acrylic resin

For fine corrections in acrylic resin, the sealing of micro-gaps and much more.

picopreci plus is a light-curing precision resin for use in the dental practice and laboratory.

The gel is ideal for fine adjustments in acrylic, sealing of micro-gaps, quick reattachment of teeth on dentures, production of denture saddles, implant edge design, precise work directly in the mouth, e.g. for the closing of screw channels, and many other purposes. picopreci plus is easy to process, has high stability and excellent modelling properties.

picopreci plus is available in pink, pink opaque (34), tooth colour and transparent.

### Further advantages:

- High stability
- Very economical dosage
- Good moulding properties

### Areas of application:

- Reparation and supplementations and correction of adjusted bite splints
- Positioning and fixing of functional elements of ortho-units, x-ray and drilling templates
- Fabrication of individual fixing and implant keys
- Fine corrections on acrylics (e.g. papillae)
- Design of papillae
- Repair of transitional prostheses
- Enables precise work in the oral cavity, e.g. sealing of screw channels

Item no.	Item	Contents
1820 0001	picopreci plus Set	3 x 3 g picopreci plus pink in light-protected syringes incl. 12 application tips and 5 ml picopreci plus Bonder

Item no.	Item	Contents
1820 0002	picopreci plus Set	3 x 3 g picopreci plus pink opaque (34) in light-protected syringes incl. 12 application tips and 5 ml picopreci plus Bonder
1820 0006	picopreci plus Set	3 x 3 g picopreci plus tooth colour in light-protected syringes incl. 12 application tips and 5 ml picopreci plus Bonder
1820 0007	picopreci plus Set	3 x 3 g picopreci plus transparent in light-protected syringes incl. 12 application tips and 5 ml picopreci plus Bonder
1820 0008	picopreci plus refill pack	5 x 3 g picopreci plus pink in light-protected syringes incl. 20 application tips
1820 0009	picopreci plus refill pack	5 x 3 g picopreci plus transparent in light-protected syringes incl. 20 application tips
1820 0005	picopreci plus Bonder	5 ml bottle
1809 0012	picopreci plus Dosierspitzen	refill pack: 50 application tips





### picobello form

#### Light-curing modelling resin

picobello is a light-curing modelling resin for making precise attachments, telescope and conical works, super-structures, implant abutments, etc.

picobello is available as a gel and now also as picobello form modelling material for swifter work with the application of larger material volumes, e.g. when making fixation and transfer keys.

When applying with the syringe or modelling by hand, you always get the right consistency for every application.

picobello gel and picobello form are compatible and can therefore be used well together.

Item no.	Item	Contents
1820 0025	picobello form	25 g can



### picobloc

#### Light-curing block-out material

picobloc gel:  
light-curing material for all types of block-outs on the model, particularly on orthodontics

picobloc wax:  
light-curing blocking-out material with high thermal stability

Item no.	Item	Contents
1809 0003	picobloc gel rot (red)	3 x 3 g syringes incl. 12 application tips
1809 0012	picobloc Dosierspitzen	refill pack: 50 application tips



### pico-light LED

#### LED polymerisation lamp

The pico-light LED is a leading-edge, high-performance light-curing device for curing light-curing precision materials in the laboratory

Perfect curing in a matter of seconds, designed for use with all light-curing precision materials from picodent®, simple operation and extremely high-performance (allows more than three hours of continuous operation).

Item no.	Item	Contents
1821 0000	pico-light LED	1 pico-light LED handpiece with light guide, 1 rechargeable battery, 1 charger, 1 safety goggles, 1 eye protection, 1 LED-Ständer (Acryl)
1821 0052	light guide for pico-light LED	1 piece
1821 0053	Rechargeable battery for pico-light LED 4.200 mAh	1 piece
1821 0054	charger for pico-light LED	1 piece
1821 0055	Safety goggles	1 piece
1821 0056	Eye protection with rubber insert ring	1 piece
1821 0010	LED-Halter (Metal)	1 piece LED-holder (Metal)
1821 0020	LED-Ständer (Acryl)	1 piece LED-stand (Acryl)



## pico-tray

Light curing tray material

To make individual impression-trays and bite templates. Fits upper and lower jaws. Available in 3 different colours.

Item no.	Item	Colour	Contents
1810 0001	pico-tray	pink	50 pieces
1810 0002	pico-tray	blue	50 pieces
1810 0003	pico-tray	white	50 pieces