



DISCOVER

NEW DIMENSIONS

Discover the new dental arch stone with fascinating properties. For extremely fine processing - with maximum hardness as a result.

pico-rock[®] select 

- ✓ Feel the difference
- ✓ Creamy consistency
- ✓ Fast and high final hardness
- ✓ Brilliant surface
- ✓ Fresh modern colours
- ✓ Scannable
- ✓ Low expansion

picodent[®] Dental-Produktions- und Vertriebs-GmbH

Lüdenscheider Str. 24-26
51688 Wipperfürth

+49 2267 6580 0
+49 2267 6580 30

✉ export@picodent.de
🌐 www.picodent.de

pico-rock® select

Experience a new dimension in dental arch stone.

50 independent laboratories accompanied pico-rock® select nationwide in Germany during the development phase and came to the conclusion, pico-rock® select feels notably different! The unique quality is the result of a newly developed manufacturing process. The creamy consistency, the perfectly matched grain structure and the particularly high fineness will inspire users.

pico-rock® select not only offers convincing processing but also achieves a high final hardness after a very short time. It can already be demoulded after approx. 35-40 min. and has a brilliant finish.

Select your favourite from the four fresh colour shades and experience the change!



Technical data

Processing time:	6 - 7 min.
Stirring time under vacuum:	45 sec.
Demoulding after:	approx. 35-40 min.
Hardness after 24h:	approx. 330 MPa (N/mm ²)
Compressive strength after 24h:	approx. 70 MPa (N/mm ²)
Expansion after 2h:	< 0.08 %

ORDER NOW AT SPECIAL PRICES

5 KG CLIMATE BAGS



1 climate bag à 5 kg
pico-rock® select powder

20 KG CARTON



4 climate bags à 5 kg
pico-rock® select powder

CHOISE OF COLOURS



pearl



toffee



sunshine



moon

Discover pico-rock® select now,
the new dental arch stone with
fascinating properties.

We look forward to your contacting us.

TELEPHONE
+49 2267 6580 0

FAX
+49 2267 6580 30

picodent® Dental-Produktions- und Vertriebs-GmbH

Lüdenscheider Str. 24-26
51688 Wipperfurth

+49 2267 6580 0
+49 2267 6580 30

export@picodent.de
www.picodent.de

